

Mouse CHEM (Chemerin) Sandwich ELISA Kit Ready-To-Use (STJE0007020)

STJE0007020

GENERAL INFORMATION

Product Type	Kits
Short Description	This CHEM Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	78.1-5000pg/mL
Sensitivity	26pg/mL
Detection Limit	78.1-5000pg/mL
Assay Time	3 hrs
Storage	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Gene	Rarres2
Symbol	
Uniprot ID	RARR2_MOUSE