

Mouse PF4 (Platelet Factor 4) Sandwich ELISA Kit Ready-To-Use (STJE0007025)

STJE0007025

GENERAL INFORMATION

Product Type	Kits
Short Description	This PF4 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	78.1-5000pg/mL
Sensitivity	27pg/mL
Detection Limit	78.1-5000pg/mL
Assay Time	3 hrs
Storage	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Instruction	table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, platelet-poor plasma, cell culture supernates or other biological fluids.
Gene	Pf4
Symbol	
Uniprot ID	PLF4_MOUSE