

Mouse RBP4 (Retinol Binding Protein 4, Plasma) Sandwich ELISA Kit Ready-To-Use (STJE0007028)

STJE0007028

GENERAL INFORMATION

Product Type	Kits
Short Description	This RBP4 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	3.12-200ng/mL
Sensitivity	1.15ng/mL
Detection Limit	3.12-200ng/mL
Assay Time	3 hrs
Storage Instruction	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	tissue homogenates or other biological fluids.
Gene	Rbp4
Symbol	
Uniprot ID	RET4_MOUSE