

Human THBS2 (Thrombospondin 2) Sandwich ELISA Kit Ready-To-Use (STJE0007030)

STJE0007030

GENERAL INFORMATION

Product Type	Kits
Short Description	This THBS2 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	62.5-4000pg/mL
Sensitivity	25.6pg/mL
Detection Limit	62.5-4000pg/mL
Assay Time	3 hrs
Storage	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma or other biological fluids.
Gene	THBS2
Symbol	
Uniprot ID	TSP2_HUMAN