

## Rat CART (Cocaine And Amphetamine Regulated Transcript) Sandwich ELISA Kit (STJE0007682)

STJE0007682

### GENERAL INFORMATION


|                          |                                                                                                                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Kits                                                                                                                                                                                     |
| <b>Short Description</b> | This CART Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA                                                                                                                                                                                    |
| <b>Reactivity</b>        | Rat                                                                                                                                                                                      |

### PRODUCT PROPERTIES

|                        |                                                                                                                                                                  |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Detection</b>       | 15.6-1000pg/mL                                                                                                                                                   |
| <b>Sensitivity</b>     | 5.8pg/mL                                                                                                                                                         |
| <b>Detection Limit</b> | 15.6-1000pg/mL                                                                                                                                                   |
| <b>Assay Time</b>      | 4,5 hrs                                                                                                                                                          |
| <b>Storage</b>         | Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| <b>Note</b>            | FOR RESEARCH USE ONLY (RUO).                                                                                                                                     |

### TARGET INFORMATION

|                     |                                                               |
|---------------------|---------------------------------------------------------------|
| <b>Sample Types</b> | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Gene</b>         | <a href="#">Cartpt</a>                                        |
| <b>Symbol</b>       |                                                               |
| <b>Uniprot ID</b>   | <a href="#">CART_RAT</a>                                      |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081