

Human CCBE1 (Collagen and calcium-binding EGF domain-containing protein 1) Sandwich ELISA Kit (STJE0007737)

STJE0007737

GENERAL INFORMATION


| | |
|--------------------------|---|
| Product Type | Kits |
| Short Description | This CCBE1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Detection | 31.2-2000pg/mL |
| Sensitivity | 11.73pg/mL |
| Detection Limit | 31.2-2000pg/mL |
| Assay Time | 4.5 hrs |
| Storage Instruction | Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|--|
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Gene | CCBE1 |
| Symbol | |
| Uniprot ID | CCBE1_HUMAN |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081