

## Mouse H4 (Histone H4) Sandwich ELISA Kit (STJE0009379)

STJE0009379

### GENERAL INFORMATION


<b>Product Type</b>	Kits
<b>Short Description</b>	This H4 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
<b>Applications</b>	ELISA
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Detection</b>	0.156-10ng/mL
<b>Sensitivity</b>	0.049ng/mL
<b>Detection Limit</b>	0.156-10ng/mL
<b>Assay Time</b>	4.5 hrs
<b>Storage</b>	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
<b>Instruction</b>	
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Gene</b>	<a href="#">H4c1</a> , <a href="#">H4c2</a> , <a href="#">H4c3</a> , <a href="#">H4c4</a> , <a href="#">H4c4</a> , <a href="#">H4c6</a> , <a href="#">H4c8</a> , <a href="#">H4c9</a> , <a href="#">H4c11</a> , <a href="#">H4c12</a> , <a href="#">Hist1h4m</a> , <a href="#">H4c14</a> , <a href="#">H4f16</a>
<b>Symbol</b>	
<b>Uniprot ID</b>	<a href="#">H4_MOUSE</a>



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081