

Rat TINAGL1 (Rat Tubulointerstitial Nephritis Antigen-Like) Sandwich ELISA Kit (STJE0012611)

STJE0012611

GENERAL INFORMATION

Product Type	Kits
Short Description	This TINAGL1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection	78.1-5000pg/mL
Sensitivity	32.6pg/mL
Detection Limit	78.1-5000pg/mL
Assay Time	4.5 hrs
Storage	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Gene	Tinagl1
Symbol	
Uniprot ID	TINAL_RAT