

## Human ASTL (Astacin Like Metallo Endopeptidase) Sandwich ELISA Kit Ready-To-Use (STJE0013730)

STJE0013730

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Kits   |
| <b>Short Description</b> | This ASTL Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA  |
| <b>Reactivity</b>        | Human  |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Detection</b>           | 0.156-10ng/mL  |
| <b>Sensitivity</b>         | 0.060ng/mL   |
| <b>Detection Limit</b>     | 0.156-10ng/mL  |
| <b>Assay Time</b>          | 3 hrs  |
| <b>Storage Instruction</b> | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| <b>Note</b>                | FOR RESEARCH USE ONLY (RUO).   |

### TARGET INFORMATION

|                     |  |
|---------------------|--|
| <b>Sample Types</b> | tissue homogenates or other biological fluids. |
| <b>Gene</b>         | <a href="#">ASTL</a>                           |
| <b>Symbol</b>       |  |
| <b>Uniprot ID</b>   | <a href="#">ASTL_HUMAN</a>                     |