

Human BRCA1 (Breast Cancer Susceptibility Protein 1) Sandwich ELISA Kit Ready-To-Use (STJE0013921)

STJE0013921

GENERAL INFORMATION

Product Type	Kits
Short Description	This BRCA1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection Sensitivity	31.2-2000pg/mL
Detection Limit	15.6pg/mL
Assay Time	31.2-2000pg/mL
Storage Instruction	3 hrs
Note	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	tissue homogenates, cell lysates or other biological fluids.
Gene Symbol	BRCA1
Uniprot ID	BRCA1_HUMAN