

Universal BSA (Bovine Serum Albumin) Competitive ELISA Kit Ready-To-Use (STJE0013926)

STJE0013926

GENERAL INFORMATION

Product Type	Kits
Short Description	This BSA Competitive ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in universal cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Universal

PRODUCT PROPERTIES

Detection	0.122-500 Mu g/mL
Sensitivity	0.05 Mu g/mL
Detection Limit	0.122-500 Mu g/mL
Assay Time	2.5 hrs
Storage Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Gene Symbol	
Uniprot ID	