

Rat C4BPα (C4 Binding Protein Alpha) Sandwich ELISA Kit Ready-To-Use (STJE0013970)

STJE0013970

GENERAL INFORMATION

Product Type	Kits
Short Description	This C4BPα Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection	62.5-4000pg/mL
Sensitivity	25.6pg/mL
Detection Limit	62.5-4000pg/mL
Assay Time	3 hrs
Storage	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Gene	C4bpa
Symbol	
Uniprot ID	C4BPA_RAT