

Rat cADPRH (Cyclic ADP Ribose Hydrolase) Sandwich ELISA Kit Ready-To-Use (STJE0014020)

STJE0014020

GENERAL INFORMATION

Product Type	Kits
Short Description	This cADPRH Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection Sensitivity	3.12-200ng/mL
Detection Limit	1.24ng/mL
Assay Time	3.12-200ng/mL
Storage Instruction	3 hrs
Note	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Gene Symbol	Cd38
Uniprot ID	CD38_RAT