

## Human CALPRO (Calprotectin) Sandwich ELISA Kit Ready-To-Use (STJE0014027)

STJE0014027

### GENERAL INFORMATION

|                     |  |
|---------------------|--|
| <b>Product Type</b> | Kits   |
| <b>Description</b>  | This CALPRO Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b> | ELISA  |
| <b>Reactivity</b>   | Human  |

### PRODUCT PROPERTIES

|                        |  |
|------------------------|--|
| <b>Detection</b>       | 31.2-2000pg/mL   |
| <b>Sensitivity</b>     | 13.12pg/mL   |
| <b>Detection Limit</b> | 31.2-2000pg/mL   |
| <b>Assay Time</b>      | 3 hrs  |
| <b>Storage</b>         | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| <b>Note</b>            | FOR RESEARCH USE ONLY (RUO).   |

### TARGET INFORMATION

|                     |  |
|---------------------|--|
| <b>Sample Types</b> | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Gene</b>         | <a href="#">S100A9</a>   |
| <b>Symbol</b>       | <a href="#">S100A8</a>   |
| <b>Uniprot ID</b>   | <a href="#">S10A9_HUMAN</a><br><a href="#">S10A8_HUMAN</a>   |