

Guinea Pig CCK8 (pig Cholecystokinin 8, Octapeptide) Competitive ELISA Kit Ready-To-Use (STJE0014121)

STJE0014121

GENERAL INFORMATION

Product Type	Kits
Short Description	This CCK8 Competitive ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in guinea pig cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Guinea Pig

PRODUCT PROPERTIES

Detection	15.6-1000pg/mL
Sensitivity	5.28pg/mL
Detection Limit	15.6-1000pg/mL
Assay Time	2.5 hrs
Storage Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates, cell culture supernates or other biological fluids.
Gene Symbol	
Uniprot ID	