

Mouse CGRP (Calcitonin Gene Related Peptide) Competitive ELISA Kit Ready-To-Use (STJE0014286)

STJE0014286

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Kits |
| Short Description | This CGRP Competitive ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Detection | 15.6-1000pg/mL |
| Sensitivity | 4.22pg/mL |
| Detection Limit | 15.6-1000pg/mL |
| Assay Time | 2.5 hrs |
| Storage Instruction | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|--|
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Gene Symbol | Calcb |
| Uniprot ID | CALCB_MOUSE |