

Human DEFA3 (Defensin Alpha 3, Neutrophil Specific) Sandwich ELISA Kit Ready-To-Use (STJE0014772)

STJE0014772

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Kits |
| Short Description | This DEFA3 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Detection | 46.8-3000pg/mL |
| Sensitivity | 16.47pg/mL |
| Detection Limit | 46.8-3000pg/mL |
| Assay Time | 3 hrs |
| Storage Instruction | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|--|
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Gene | DEFA3 |
| Symbol | |
| Uniprot ID | DEF3_HUMAN |