

Rat DEFb2 (Defensin Beta 2) Sandwich ELISA Kit Ready-To-Use (STJE0014792)

STJE0014792

GENERAL INFORMATION

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type | Kits |
| Short Description | This DEFb2 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Detection | 0.156-10ng/mL |
| Sensitivity | 0.061ng/mL |
| Detection Limit | 0.156-10ng/mL |
| Assay Time | 3 hrs |
| Storage | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|---------------------------------------------------------------|
| Sample Types | serum, plasma, tissue homogenates or other biological fluids. |
| Gene | Defb4 |
| Symbol | |
| Uniprot ID | DEFB4_RAT |