

Human FcγR3A (Fc Fragment Of IgG Low Affinity IIIa Receptor) Sandwich ELISA Kit Ready-To-Use (STJE0015206)

STJE0015206

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This FcγR3A Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|------------------------|--|
| Detection | 78.1-5000ng/mL |
| Sensitivity | 34ng/mL |
| Detection Limit | 78.1-5000ng/mL |
| Assay Time | 3 hrs |
| Storage | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | serum, plasma, tissue homogenates, cell lysates or other biological fluids. |
| Gene | FCGR3A |
| Symbol | |
| Uniprot ID | FCG3A_HUMAN |