

Mouse GFER (Growth Factor, Augmenter Of Liver Regeneration) Sandwich ELISA Kit Ready-To-Use (STJE0015517)

STJE0015517

GENERAL INFORMATION

Product Type	Kits
Short Description	This GFER Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection Sensitivity	78.1-5000pg/mL
Detection Limit	29pg/mL
Assay Time	78.1-5000pg/mL
Storage	3 hrs
Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma or other biological fluids.
Gene Symbol	Gfer
Uniprot ID	ALR_MOUSE