

## Rat GSTt2 (Glutathione S Transferase Theta 2) Sandwich ELISA Kit Ready-To-Use (STJE0015732)

STJE0015732

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This GSTt2 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
<b>Applications</b>	ELISA
<b>Reactivity</b>	Rat

### PRODUCT PROPERTIES

<b>Detection Sensitivity</b>	1.56-100U/L
<b>Detection Limit</b>	1U/L
<b>Assay Time</b>	1.56-100U/L
<b>Storage</b>	3 hrs
<b>Instruction</b>	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Gene</b>	<a href="#">Gstt2</a>
<b>Symbol</b>	
<b>Uniprot ID</b>	<a href="#">GSTT2_RAT</a>