

Mouse GUCA2A (Guanylate Cyclase Activator 2A) Sandwich ELISA Kit Ready-To-Use (STJE0015740)

STJE0015740

GENERAL INFORMATION

Product Type	Kits
Short Description	This GUCA2A Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection Sensitivity	31.2-2000fmol/mL
Detection Limit	8.45fmol/mL
Assay Time	31.2-2000fmol/mL
Storage	3 hrs
Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Gene	Guca2a
Symbol	
Uniprot ID	GUC2A_MOUSE