

## Mouse H4 (Histone H4) Sandwich ELISA Kit Ready-To-Use (STJE0015760)

STJE0015760

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This H4 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
<b>Applications</b>	ELISA
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Detection</b>	0.156-10ng/mL
<b>Sensitivity</b>	0.049ng/mL
<b>Detection Limit</b>	0.156-10ng/mL
<b>Assay Time</b>	3 hrs
<b>Storage</b>	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Gene</b>	<a href="#">H4c1.H4c2.H4c3.H4c4.H4c6.H4c8.H4c9.H4c11.H4c12.Hist1h4m.H4c14.H4f16</a>
<b>Symbol</b>	
<b>Uniprot ID</b>	<a href="#">H4_MOUSE</a>