

Rabbit HCII (Heparin Cofactor II) Sandwich ELISA Kit Ready-To-Use (STJE0015790)

STJE0015790

GENERAL INFORMATION

Product Type	Kits
Description	This HCII Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rabbit cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Rabbit

PRODUCT PROPERTIES

Detection	3.12-200U/L
Sensitivity	1.03U/L
Detection Limit	3.12-200U/L
Assay Time	3 hrs
Storage	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	plasma.
Gene	
Symbol	
Uniprot ID	