

## Mouse HMGCR (3-Hydroxy-3-Methylglutaryl Coenzyme A Reductase) Sandwich ELISA Kit Ready-To-Use (STJE0015860)

STJE0015860

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Kits  |
| <b>Short Description</b> | This HMGCR Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA   |
| <b>Reactivity</b>        | Mouse   |

### PRODUCT PROPERTIES

|                        |  |
|------------------------|--|
| <b>Detection</b>       | 0.781-50ng/mL  |
| <b>Sensitivity</b>     | 0.31ng/mL  |
| <b>Detection Limit</b> | 0.781-50ng/mL  |
| <b>Assay Time</b>      | 3 hrs  |
| <b>Storage</b>         | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| <b>Note</b>            | FOR RESEARCH USE ONLY (RUO).   |

### TARGET INFORMATION

|                     |  |
|---------------------|--|
| <b>Sample Types</b> | tissue homogenates, cell lysates or other biological fluids. |
| <b>Gene</b>         | <a href="#">Hmgcr</a>  |
| <b>Symbol</b>       |  |
| <b>Uniprot ID</b>   | <a href="#">HMDH_MOUSE</a>                                   |