

Mouse IL10Ra (Interleukin 10 Receptor Alpha) Sandwich ELISA Kit Ready-To-Use (STJE0016134)

STJE0016134

GENERAL INFORMATION

Product Type	Kits
Short Description	This IL10Ra Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	78.1-5000pg/mL
Sensitivity	29pg/mL
Detection Limit	78.1-5000pg/mL
Assay Time	3 hrs
Storage Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	tissue homogenates, cell lysates or other biological fluids.
Gene	Il10ra
Symbol	
Uniprot ID	I10R1_MOUSE