

## Mouse ITIH4 (Inter Alpha-Globulin Inhibitor H4) Sandwich ELISA Kit Ready-To-Use (STJE0016486)

STJE0016486

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This ITIH4 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
<b>Applications</b>	ELISA
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Detection</b>	78.1-5000ng/mL
<b>Sensitivity</b>	33ng/mL
<b>Detection Limit</b>	78.1-5000ng/mL
<b>Assay Time</b>	3 hrs
<b>Storage</b>	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Gene</b>	<a href="#">Itih4</a>
<b>Symbol</b>	
<b>Uniprot ID</b>	<a href="#">ITIH4_MOUSE</a>