

## Human JAM2 (Junctional Adhesion Molecule 2) Sandwich ELISA Kit Ready-To-Use (STJE0016506)

STJE0016506

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This JAM2 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
<b>Applications</b>	ELISA
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Detection Sensitivity</b>	31.2-2000pg/mL
<b>Detection Limit</b>	12.6pg/mL
<b>Assay Time</b>	31.2-2000pg/mL
<b>Storage Instruction</b>	3 hrs
<b>Note</b>	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Gene Symbol</b>	JAM2
<b>Uniprot ID</b>	JAM2_HUMAN