

Human MAPRE1 (Microtubule Associated Protein RP/EB Family, Member 1) Sandwich ELISA Kit Ready-To-Use (STJE0016873)

STJE0016873

GENERAL INFORMATION

Product Type	Kits
Short Description	This MAPRE1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	0.156-10ng/mL
Sensitivity	0.036ng/mL
Detection Limit	0.156-10ng/mL
Assay Time	3 hrs
Storage Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	tissue homogenates, cell lysates or other biological fluids.
Gene	MAPRE1
Symbol	
Uniprot ID	MARE1_HUMAN