

Human MMRN2 (Multimerin 2) Sandwich ELISA Kit Ready-To-Use (STJE0017087)

STJE0017087

GENERAL INFORMATION

Product Type	Kits
Short Description	This MMRN2 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	78.1-5000pg/mL
Sensitivity	33pg/mL
Detection Limit	78.1-5000pg/mL
Assay Time	3 hrs
Storage	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Gene	MMRN2
Symbol	
Uniprot ID	MMRN2_HUMAN