

Rat MYH7 (Myosin Heavy Chain 7, Cardiac Muscle, Beta) Sandwich ELISA Kit Ready-To-Use (STJE0017215)

STJE0017215

GENERAL INFORMATION

Product Type	Kits
Short Description	This MYH7 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection	125-8000pg/mL
Sensitivity	43pg/mL
Detection Limit	125-8000pg/mL
Assay Time	3 hrs
Storage	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Gene	Myh7
Symbol	
Uniprot ID	MYH7_RAT