

Human NIS (Sodium Iodide Symporter) Sandwich ELISA Kit Ready-To-Use (STJE0017361)

STJE0017361

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Kits |
| Short Description | This NIS Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|------------------------|--|
| Detection | 0.781-50ng/mL |
| Sensitivity | 0.24ng/mL |
| Detection Limit | 0.781-50ng/mL |
| Assay Time | 3 hrs |
| Storage | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | serum, plasma, tissue homogenates, cell lysates or other biological fluids. |
| Gene | SLC5A5 |
| Symbol | |
| Uniprot ID | SC5A5_HUMAN |