

Rabbit PAI1 (Plasminogen Activator Inhibitor 1) Sandwich ELISA Kit Ready-To-Use (STJE0017643)

STJE0017643

GENERAL INFORMATION

Product Type	Kits
Short Description	This PAI1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rabbit cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Rabbit

PRODUCT PROPERTIES

Detection	1.56-100pg/mL
Sensitivity	0.68pg/mL
Detection Limit	1.56-100pg/mL
Assay Time	3 hrs
Storage Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	plasma.
Gene Symbol	
Uniprot ID	