

Mouse PDE3B (Phosphodiesterase 3B, cGMP Inhibited) Sandwich ELISA Kit Ready-To-Use (STJE0017731)

STJE0017731

GENERAL INFORMATION

Product Type	Kits
Short Description	This PDE3B Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	31.2-2000pg/mL
Sensitivity	18.75pg/mL
Detection Limit	31.2-2000pg/mL
Assay Time	3 hrs
Storage Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Gene	Pde3b
Symbol	
Uniprot ID	PDE3B_MOUSE