

Rat PDK2 (Pyruvate Dehydrogenase Kinase Isozyme 2) Sandwich ELISA Kit Ready-To-Use (STJE0017767)

STJE0017767

GENERAL INFORMATION

Product Type	Kits
Short Description	This PDK2 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection Sensitivity	78.1-5000pg/mL
Detection Limit	35.6pg/mL
Assay Time	78.1-5000pg/mL
Storage	3 hrs
Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Gene	Pdk2
Symbol	
Uniprot ID	PDK2_RAT