

## Rat PGC (Pepsinogen C) Sandwich ELISA Kit Ready-To-Use (STJE0017813)

STJE0017813

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This PGC Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
<b>Applications</b>	ELISA
<b>Reactivity</b>	Rat

### PRODUCT PROPERTIES

<b>Detection</b>	4.68-300ng/mL
<b>Sensitivity</b>	1.35ng/mL
<b>Detection Limit</b>	4.68-300ng/mL
<b>Assay Time</b>	3 hrs
<b>Storage</b>	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Gene</b>	<a href="#">Pgc</a>
<b>Symbol</b>	
<b>Uniprot ID</b>	<a href="#">PEPC_RAT</a>