

Mouse PP2Aa (Serine/Threonine-Protein Phosphatase 2A 56 kDa Regulatory Subunit Alpha Isoform) Sandwich ELISA Kit Ready-To-Use (STJE0017988)

STJE0017988

GENERAL INFORMATION

Product Type	Kits
Short Description	This PP2Aa Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	78.1-5000pg/mL
Sensitivity	34.6pg/mL
Detection Limit	78.1-5000pg/mL
Assay Time	3 hrs
Storage	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Gene	Ppp2r1a
Symbol	
Uniprot ID	2AAA_MOUSE