

## Rat PPAR $\gamma$ C1a (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) Sandwich ELISA Kit Ready-To-Use (STJE0018001)

STJE0018001

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This PPAR $\gamma$ C1a Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
<b>Applications</b>	ELISA
<b>Reactivity</b>	Rat

### PRODUCT PROPERTIES

<b>Detection</b>	0.156-10ng/mL
<b>Sensitivity</b>	0.061ng/mL
<b>Detection Limit</b>	0.156-10ng/mL
<b>Assay Time</b>	3 hrs
<b>Storage</b>	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

<b>Sample Types</b>	tissue homogenates, cell lysates or other biological fluids.
<b>Gene</b>	<a href="#">Ppargc1a</a>
<b>Symbol</b>	
<b>Uniprot ID</b>	<a href="#">PRGC1_RAT</a>