

Human PPP3R1 (Protein Phosphatase 3, Regulatory Subunit 1) Sandwich ELISA Kit Ready-To-Use (STJE0018014)

STJE0018014

GENERAL INFORMATION

Product Type	Kits
Short Description	This PPP3R1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	31.2-2000pg/mL
Sensitivity	15.5pg/mL
Detection Limit	31.2-2000pg/mL
Assay Time	3 hrs
Storage Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	tissue homogenates or other biological fluids.
Gene	PPP3R1
Symbol	
Uniprot ID	CANB1_HUMAN