

## Human PVRL1 (Poliovirus Receptor Related Protein 1) Sandwich ELISA Kit Ready-To-Use (STJE0018165)

STJE0018165

### GENERAL INFORMATION

|                          |                                                                                                                                                                                                           |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Kits                                                                                                                                                                                                      |
| <b>Short Description</b> | This PVRL1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA                                                                                                                                                                                                     |
| <b>Reactivity</b>        | Human                                                                                                                                                                                                     |

### PRODUCT PROPERTIES

|                        |                                                                                                                                                                                                                                                |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Detection</b>       | 0.156-10ng/mL                                                                                                                                                                                                                                  |
| <b>Sensitivity</b>     | 0.063ng/mL                                                                                                                                                                                                                                     |
| <b>Detection Limit</b> | 0.156-10ng/mL                                                                                                                                                                                                                                  |
| <b>Assay Time</b>      | 3 hrs                                                                                                                                                                                                                                          |
| <b>Storage</b>         | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| <b>Note</b>            | FOR RESEARCH USE ONLY (RUO).                                                                                                                                                                                                                   |

### TARGET INFORMATION

|                     |                                                              |
|---------------------|--------------------------------------------------------------|
| <b>Sample Types</b> | tissue homogenates, cell lysates or other biological fluids. |
| <b>Gene</b>         | <a href="#">NECTIN1</a>                                      |
| <b>Symbol</b>       |                                                              |
| <b>Uniprot ID</b>   | <a href="#">NECT1_HUMAN</a>                                  |