

Human RIPK3 (Receptor Interacting Serine Threonine Kinase 3) Sandwich ELISA Kit Ready-To-Use (STJE0018280)

STJE0018280

GENERAL INFORMATION

Product Type	Kits
Short Description	This RIPK3 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection Sensitivity	31.2-2000pg/mL
Detection Limit	11.7pg/mL
Assay Time	31.2-2000pg/mL
Storage	3 hrs
Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	tissue homogenates, cell lysates or other biological fluids.
Gene Symbol	RIPK3
Uniprot ID	RIPK3_HUMAN