

Universal rT3 (Reverse Triiodothyronine) Competitive ELISA Kit Ready-To-Use (STJE0018339)

STJE0018339

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This rT3 Competitive ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in universal cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Universal |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Detection | 93.7-6000pg/mL |
| Sensitivity | 24.32pg/mL |
| Detection Limit | 93.7-6000pg/mL |
| Assay Time | 2.5 hrs |
| Storage Instruction | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | serum, plasma or other biological fluids. |
| Gene Symbol | |
| Uniprot ID | |