

Human SH3L3 (SH3 Domain-Binding Glutamic Acid-Rich-Like Protein 3 (SH3L3)) Sandwich ELISA Kit Ready-To-Use (STJE0018536)

STJE0018536

GENERAL INFORMATION

Product Type	Kits
Short Description	This SH3L3 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	78.1-5000pg/mL
Sensitivity	33.4pg/mL
Detection Limit	78.1-5000pg/mL
Assay Time	3 hrs
Storage	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Gene	SH3BGRL3
Symbol	
Uniprot ID	SH3L3_HUMAN