

Guinea Pig TGFb1 (pig Transforming Growth Factor Beta 1) Sandwich ELISA Kit Ready-To-Use (STJE0018898)

STJE0018898

GENERAL INFORMATION

Product Type	Kits
Short Description	This TGFb1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in guinea pig cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Guinea Pig

PRODUCT PROPERTIES

Detection	31.2-2000pg/mL
Sensitivity	7.6pg/mL
Detection Limit	31.2-2000pg/mL
Assay Time	3 hrs
Storage	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Gene	
Symbol	
Uniprot ID	