

Equine TIMP1 (Tissue Inhibitors Of Metalloproteinase 1) Sandwich ELISA Kit Ready-To-Use (STJE0018970)

STJE0018970

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This TIMP1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in equine cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Equine |

PRODUCT PROPERTIES

| | |
|------------------------------|--|
| Detection Sensitivity | 62.5-4000pg/mL |
| Detection Limit | 25.3pg/mL |
| Assay Time | 62.5-4000pg/mL |
| Storage | 3 hrs |
| Instruction | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | serum, plasma, tissue homogenates or other biological fluids. |
| Gene Symbol | |
| Uniprot ID | |