

## Rat TINAGL1 (Rat Tubulointerstitial Nephritis Antigen-Like) Sandwich ELISA Kit Ready-To-Use (STJE0018999)

STJE0018999

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Kits  |
| <b>Short Description</b> | This TINAGL1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA   |
| <b>Reactivity</b>        | Rat   |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Detection</b>           | 78.1-5000pg/mL   |
| <b>Sensitivity</b>         | 32.6pg/mL  |
| <b>Detection Limit</b>     | 78.1-5000pg/mL   |
| <b>Assay Time</b>          | 3 hrs  |
| <b>Storage Instruction</b> | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| <b>Note</b>                | FOR RESEARCH USE ONLY (RUO).   |

### TARGET INFORMATION

|                     |   |
|---------------------|---|
| <b>Sample Types</b> | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Gene</b>         | <a href="#">Tinagl1</a>                                       |
| <b>Symbol</b>       |   |
| <b>Uniprot ID</b>   | <a href="#">TINAL_RAT</a>                                     |