

Human TMSb10 (Thymosin Beta 10) Competitive ELISA Kit Ready-To-Use (STJE0019048)

STJE0019048

GENERAL INFORMATION

Product Type	Kits
Short Description	This TMSb10 Competitive ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	7.81-500ng/mL
Sensitivity	3.52ng/mL
Detection Limit	7.81-500ng/mL
Assay Time	2.5 hrs
Storage	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Gene	TMSB10
Symbol	
Uniprot ID	TYB10_HUMAN