

Mouse TNFSF13 (Tumor Necrosis Factor Ligand Superfamily, Member 13) Sandwich ELISA Kit Ready-To-Use (STJE0019086)

STJE0019086

GENERAL INFORMATION

Product Type	Kits
Short Description	This TNFSF13 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	46.8-3000pg/mL
Sensitivity	22.7pg/mL
Detection Limit	46.8-3000pg/mL
Assay Time	3 hrs
Storage Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma or other biological fluids.
Gene	Tnfsf13
Symbol	
Uniprot ID	TNF13_MOUSE