

Rat TNNI1 (Troponin I Type 1, Slow Skeletal) Sandwich ELISA Kit Ready-To-Use (STJE0019098)

STJE0019098

GENERAL INFORMATION

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type | Kits |
| Short Description | This TNNI1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in rat cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Detection | 31.2-2000pg/mL |
| Sensitivity | 10.7pg/mL |
| Detection Limit | 31.2-2000pg/mL |
| Assay Time | 3 hrs |
| Storage | The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit. |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|-------------------------------------------|
| Sample Types | serum, plasma or other biological fluids. |
| Gene | Tnni1 |
| Symbol | |
| Uniprot ID | TNNI1_RAT |