

Mouse TRPM6 (Transient Receptor Potential Cation Channel Subfamily M, Member 6) Sandwich ELISA Kit Ready-To-Use (STJE0019208)

STJE0019208

GENERAL INFORMATION

Product Type	Kits
Short Description	This TRPM6 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection Sensitivity	62.5-4000pg/mL
Detection Limit	27.6pg/mL
Assay Time	62.5-4000pg/mL
Storage	3 hrs
Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	tissue homogenates or other biological fluids.
Gene	Trpm6
Symbol	
Uniprot ID	TRPM6_MOUSE